

ComforGel Liner

The ComforGel Liner is a unique design that sandwiches a special formulated gel between inner and outer layers of cotton fabric. The stretch fabric allows for the liner to conform to virtually any shape and provide exceptional comfort, while the nylon reinforced distal end provides for great durability. The great stretch allows for ComforGel Liners to be used for AK, BK, Arm, and even Symes or Chopart applications!

Available in three and five ply thicknesses and a variety of sizes to accommodate practically any amputee. Standard gel length is 9", however a 12" gel length is also available. X-Tra Styles have gel for the full length of the sock.

Features of the ComforGel Liner include:

- Exceptional comfort for the amputee
- Gel formulated specifically for prosthetic use
- Unique sandwich design surrounds gel with durable cotton fabric
- Standard version has 9" of gel
- Extended gel is available, up to the full length of the sock
- Standard lengths of 12" and 14"
- Longer lengths available upon request
- Two thicknesses available, easily identified by Yellow stitching for 3 ply and Green stitching for 5 ply
- Appropriate for any level of amputation
- X-Tra version is ideal for Symes and Chopart level amputations
- Very cost effective!

**Need extra comfort
with durability?**

Look no further!

**ComforGel
Liner**



ComforGel Liner

Sizes and Part Numbers

Standard Styles						Specials - Not Returnable
Type	Ply	Sock Length	Flat Width	Maximum Circumference	9" Gel Length Part #	12" Gel Length Part #
Ped.	3 ply	12"	2 ¼"	10"	CGELYP	CGELYP/FD
Ped.	3 ply	14"	2 ¼"	10"	CGELYP14	CGELYP14/12
Ped.	5 ply	12"	2 ¼"	10"	CGELGP	CGELGP/FD
Ped.	5 ply	14"	2 ¼"	10"	CGELGP14	CGELGP14/12
AK	3 ply	12"	5 ¼"	26"	CGELAY	CGELAY/FD
AK	3 ply	14"	5 ¼"	26"	CGELAY14	CGELAY14/12
AK	5 ply	12"	5 ¼"	26"	CGELAG	CGELAG/FD
AK	5 ply	14"	5 ¼"	26"	CGELAG14	CGELAG14/12
BK	3 ply	12"	4"	20"	CGELY	CGELY/FD
BK	3 ply	14"	4"	20"	CGELY14	CGELY14/12
BK	5 ply	12"	4"	20"	CGELG	CGELG/FD
BK	5 ply	14"	4"	20"	CGELG14	CGELG14/12
Arm	3 ply	12"	2 ¼"	10"	CGELYA	CGELYA/FD
Arm	3 ply	14"	2 ¼"	10"	CGELYA14	CGELYA14/12
Arm	5 ply	12"	2 ¼"	10"	CGELGA	CGELGA/FD
Arm	5 ply	14"	2 ¼"	10"	CGELGA14	CGELGA14/12

SPECIALS - NOT RETURNABLE

X-Tra Styles with Full Length Gel					
Type	Ply	Length	Flat Width	Max. Circ.	Part #
Adult	3 ply	20"	4"	20"	CGELY20FD
Adult	3 ply	22"	4"	20"	CGELY22FD
Adult	5 ply	20"	4"	20"	CGELG20FD
Adult	5 ply	22"	4"	20"	CGELG22FD
Ped.	3 ply	20"	2 ¼"	12"	CGELYP20FD
Ped.	3 ply	22"	2 ¼"	12"	CGELYP22FD
Ped.	5 ply	20"	2 ¼"	12"	CGELGP20FD
Ped.	5 ply	22"	2 ¼"	12"	CGELGP22FD

INSTRUCTIONS

CARE: Turn sock inside out. Machine wash using cool water and delicate fabric cycle. Spin dry in washer. Hang dry with open end down to drain. Do not dry in dryer.

Avoid direct sunlight. No bleach. Do not wring.

DO NOT: Wring or twist - possible damage to gel and/or fabric could result. Excessive heat could cause loss of mineral oil and degradation of gel. Bleach will cause loss of fabric strength and premature failure.

Add a "3" to the end of any ComforGel part number for a vulcanized 3S hole
Custom Sizes Available Upon Request!